


Standard Single Use Syringes



Single Use Syringes Luer TBC scale



1. The narrow silicone seal minimizes the contact area between the silicone and the solution thereby reducing solid particle contamination. The silicone seal provides easy and accurate dosage management.
2. Luer-Slip centric nozzle. TBC scale. The dual scale provides easy and precise reading irrespective of the handling of the syringe.
3. The plunger requires despite of reduced amount of silicone oil a minimum use of force due to smooth plunger action even when used single-handed. Finger grip of an ergonomic design makes the syringe easy to handle.

REF	Product Description	Packaging Units
62.1256	1 ML ONCE Luer TBC 0,01 grad Bulk Non-Sterile	6000/96000
62.1202	1 ML ONCE Luer TBC 0,01 grad	100/3200/51200
62.1656	1 ML CODAN LUER TBC 0,01 grad Bulk Non-Sterile	100/3200/51200
62.1612	1 ML CODAN LUER TBC 0,01 grad	100/3200/51200